## dispoCHEM MAX

## Rękawice ochronne kat. III

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X - Test method appears not to be suitable for the glove design/ material

**Sizes:** 7/S, 8/M, 9/L, 10/XL, 11/XXL

**Packing:** 50 pieces

AMBIDEXTROUS  
NON-STERIL

**Product description:**  
Disposable powder free glove, 100% Nitrile Latex (NBR)

Glove pH 6.0 - 8.0  
dexterity level 5

### General

100% Nitrile Latex (NBR) powder free glove, free of natural rubber, specially formulated to provide chemical protection against specified solvents, chemicals and diluted pesticides. Also this glove has good dexterity. Chlorinated for facile donning and doffing. This glove features a rolled cuff for added protection from potential chemical spills and drips and is ambidextrous. None of the materials or processes used in the manufacture of these products are known to be harmful to the wearer. The product is in compliance with the PPE Regulation (EU) 2016/425. The products meet regulation REACH 1907/2006 Annex XVII Entry number 50 requirements for PAHs. Test results are taken from palm area of the glove and protection is limited to the palm only. Glove performance quoted is based on laboratory data and may not reflect the actual duration of protection in the workplace due to the differentiation between mixtures and pure chemicals and other factors influencing the performance such as temperature, abrasion, degradation etc. The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only (except in cases where the glove is equal to or over 400 mm - where the cuff is tested also) and relates only to the chemical tested. It can be different if the chemical is used in a mixture. When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves. The duration of the test is not based on actual use time since the permeation test is an accelerated test in which the surface of the specimen is in constant contact with the test chemical. Although the duration of the exposure may be for a longer period during field application with a dilute formulation, the entire surface is not in constant contact with the test chemical. Gloves are suitable when the potential risk is relatively low. These gloves are not suitable for use with concentrated pesticide formulations and/or for scenarios where mechanical risks exist. Check the minimum overlap. If the overlap is less than approximately 50 mm between the glove and the sleeve, a glove with longer length should be used. The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen. When an indication for hand hygiene precedes a contact that also requires glove usage, hand rubbing or hand washing should be performed before donning gloves and after removing. Degradation levels indicate the change in puncture resistance after exposure to the challenge chemical.

### Glove donning

Check for any physical damages, glove condition and contamination prior to each use of glove. Ensure the selection of appropriate glove size. Tactile sensitivity and dexterity reduce with poor fitting of glove cause fatigue in hands and finger. Wrong size and poor fitting leads to poor hand protection.

### Glove doffing

Removed as soon as showing signs of wear. When contaminant is not removable or presents a potential hazard it is advisable to ease left and right gloves off alternately using the gloved hand so that the gloves are removed without the contaminant contacting bare hands.

### Obsolescence

Ideally stored in dry conditions between 10°C to 30°C in its original package. Do not exceed 40°C. When stored as recommended will not suffer change in chemical properties for up to five years from date of manufacture. Service life cannot be specified and depends on application and responsibility of user to ascertain suitability of the glove for its intended use.

### Cleaning / Maintenance

For single use only. Before usage, inspect the gloves for any defect or imperfections. Do not subjected to mechanical cleaning procedures. Gloves should be cleaned as far as possible provided no serious hazard exists, before removal from hand. When contaminant is not removable or presents a potential hazard it is advisable to ease left and right gloves off alternately using the gloved hand so that the gloves are removed without the contaminant contacting bare hands.

### Disposal

Used gloves should be disposed by incineration or as per the regulation of local authorities. **Please note:** It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation.

**Caution:** The results of the tests should help in glove selection, however it must be understood that actual conditions of use cannot be simulated and it is the responsibility of the user not the manufacturer to determine glove suitability to the intended use. Glove contains accelerators may cause allergic reactions. Discontinue use of this product in case of rashes appear on the skin and seek medical advises. Further information may be obtained from manufacturer.  
Glove is in compliance with the regulation (EC) No 1935/2004 and is suitable for handling all kinds of foodstuffs for short period of time.

Module D assessment carried out by:  
SGS Fimko Oy, Takomotti 8,  
FI-00380 Helsinki, Finland  
(Notified Body No. 05798)

The product is in compliance with  
the PPE Regulation (EU) 2016/425  
EU Type examination carried out by: SATRA Technology Europe Ltd,  
Bracetown Business Park, Clonee, Dublin, D15 YN2P, Ireland.  
(Notified Body No. 2777)

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CAT III

G1  
ISO 18889

EU Type examination carried out by:- SATRA Technology Europe Ltd, Bracetown Business Park,  
Clonee, Dublin, D15 YN2P, Ireland. (Notified Body No. 2777)

**Permeation by chemicals in accordance with EN ISO 374-1:2016+A1:2018//Type A**

Chemical	Code letter	Breakthrough time	Level	<i>The lowest Breakthrough time is used exdetermination of performance level.</i>					
Metanol	A	4 minutes	0						
n-Heptane	J	>480 minutes	6						
Sodium hydroxide 40%	K	>480 minutes	6						
Sulphuric Acid 96%	L	14 minutes	1						
Nitric Acid 65%	M	8 minutes	0						
Acetic Acid 99%	N	12 minutes	1						
Ammonium hydroxide 25%	O	39 minutes	2						
Hydrogen Peroxide 30%	P	>480 minutes	6						
Hydrofluoric Acid 40%	S	35 minutes	2						
Formaldehyde 37%	T	223 minutes	4						
<i>Performance level</i>	1	2	3	4	5	6			
<i>Breakthrough time (minutes)</i>	>10	>30	>60	>120	>240	>480			

**Micro-organisms risk in accordance with EN ISO 374-5:2016**

Protection against bacteria & fungi: **PASS**

Protection against Viruses: **PASS**

**Degradation by chemicals in accordance with EN ISO 374-4:2019**

Chemical	Code letter	Mean Degradation	Appearance of the sample after testing
Metanol	A	81,8%	Moderate swelling
n-Heptane	J	13,6%	Slight swelling
Sodium hydroxide 40%	K	5,0%	No change
Sulphuric Acid 96%	L	100%	Complete degradation
Nitric Acid 65%	M	98,5%	Sever swelling and colour change
Acetic Acid 99%	N	96,3%	Sever swelling
Ammonium hydroxide 25%	O	34,7%	Moderate swelling
Hydrogen Peroxide 30%	P	43,4%	Slight swelling
Formaldehyde 37%	T	21,7%	Slight swelling
Hydrofluoric Acid 40%	S	-	-

### Mechanical test data in accordance with EN 388:2016+A1:2018

Mechanical test data	Result	Mechanical Property / Level	1	2	3	4	5
Abrasion resistance	Level 2	<i>Abrasion resistance (Cycles)</i>	100	500	2000	8000	–
Cut resistance	Level 0	<i>Cut resistance (Index)</i>	1,2	2,5	5,0	10,0	20,0
Tear resistance	Level 0	<i>Tear resistance (Newton)</i>	10	25	50	75	–
Puncture resistance	Level 0	<i>Puncture resistance (Newton)</i>	20	60	100	150	–
TDM Cut Resistance	X	<i>Cutting force (Newton)</i>	A >2, B >5, C >10, D >15, E >22, F >30				

*X-Not tested*

*0 – Level 1 not achieved*

### Permeation by chemical in accordance with ISO 18889:2019

Chemical	Test Method	Requirement	Test Results
EC-DY (surrogate pesticide) diluted solution	ISO 19918:2017	≤10 µg/cm <sup>2</sup>	Pass

Micro-organisms risk in accordance with EN ISO 374-5:2016

Protection against bacteria & fungi: PASS

Protection against Viruses: PASS