nitrylBLUE

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3111X

ISO 18889

Product is in compliance with requirements of the PPE Regulation 2016/425

X - Test method appears not to be suitable for the glove design/material.

EU- Type Examination Certificate

 SATRA Technology Europe Limited
 Braesawn Business Park, Clonee, D15 YN2P, Dublin 15, Ireland.
 Notified Body: 2777

Precautions for use

Before using the gloves make sure that there are no mechanical damage. Damaged gloves should be replaced by new ones. Before working, you should choose the correct size of gloves. When an indication of hand hygiene precedes a contact that also requires glove usage, hand rubbing or hand washing should be performed before donning glove. Do not use in position where there is a risk of entanglement in rotating parts of the machine. Gloves meeting the requirement for resistance to puncture may not be suitable for protection against sharply pointed objects such as hypodermic needles. These gloves are tested on EC-DY 1 hours with diluted solution that provide protection against Pesticides (Chemical Hand Protection) Remove the glove immediately if contaminated by a concentrate spill. The absorbent lining material has the potential of absorbing the pesticide. This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals. It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation. When used, protective gloves may provide less resistance to the dangerous chemicals due to changes in physical properties. Movements, snagging, rubbing, degradation caused. These gloves are tested on EC-DY 1 hours with diluted solution that provide protection against Pesticides (Chemical Hand Protection) Remove the glove immediately if contaminated by a concentrate spill. The absorbent lining material has the potential of absorbing the pesticide. This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals. It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation. When used, protective gloves may provide less resistance to the dangerous chemicals due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves. GR gloves provide protection only to the palm-side of the hand for re-entry worker who is in contact with dry and partially dry pesticide residues that remain on the plant surface after pesticide application. This glove category is suitable only for re-entry activities where it has been determined that protection provided to the fingertips and palm-side of the hand is sufficient. **Doffing:** Clean off the outside of the gloves at least daily and prior to removing them. Loosen the gloves on both hands. Hold hands down to prevent the material, product, or contaminated water from dripping onto your skin or clothing. Keep the gloves down and remove the first glove only to the fingers - the cuff will remain over your palm. Grab the second glove with the first glove and remove it. The first glove should slip off. Handle gloves by the inside only - the area not contaminated. Avoid touching the outside of the glove.

Storage

Store in a cool, dry place. Do not expose to sources of ozone and open flames. Transport and storage only in supplier's packaging. Gloves should be stored at 5-25° C. When stored as recommended will not suffer change in mechanical properties for up to 5 years from the date of manufacture. Service life cannot be specified and depends on the application and responsibility to user to ascertain suitability of the glove for its intended use.

Disposal

Used gloves should be disposed by incineration or as per the regulation of local authorities. Please dispose according to local authority regulations. In case of contact with chemicals, they should be disposed of in accordance with the chemical's manufacturer's instructions.

Cleaning and care instructions

Gloves should not be left in contaminated condition if re-use is intended in which case gloves should be cleaned as far as possible provided no serious hazard exists, before removal from hand. After cleaning gloves should be thoroughly inspected. When contaminant is not removable or presents a potential hazard, it is advisable to ease left and right gloves off alternately using the gloved hand so that the gloves are removed without contaminating bare hands. Used gloves should be disposed by incineration or as per the regulation of local authorities. Gloves can be cleaned by wiping with damp cloth. This glove is not washable. The user is responsible for the protective properties of the used gloves after washing or cleaning. SEMEX does not bear any responsibility in this regard.

General

This product is manufactured under a Quality system which has been registered and meets the requirements of ISO 9001:2015. The manufacturer was examined under the system of ensuring EU Quality of production by means of monitoring EU Regulation 2016/425 Module D by SGS Fimko Oy, P.O. Box 30 (Särkinieentie 3), 00211 Helsinki, Finland (Notified Body No. 0578). The models referred to are designed to accommodate the basic safety requirements and standards laid down in EU regulation for Personal Protective Equipment and EN ISO 21420:2020, EN 388:2016+A1:2018, EN ISO 374-1:2016, EN ISO 374-5:2016 & ISO 18889:2019 respectively. None of the materials or processes used in the manufacture of these products are known to have any harmful effect on the wearer.

Please note:

Hand protection is limited to the Nitrile Rubber Coated area of the glove. The information contained herein is intended to assist the wearer in selection of Personal Protective Equipment. However, it should be understood that the glove properties stated are based on laboratory data and may not reflect the actual conditions under which gloves will be used. User is responsible for the selection of suitable gloves for the actual conditions and the length of their application. Test is taken from palm area of the glove.

LOT#14022023287

 EU Declaration of Conformity: <https://hfsafety.pl/>
Date of production: 05/2023

P.H.P.U. SEMEX Fertac, Huszno Sp. J.
 ul. Jagiellońska 101/105
 42-202 Częstochowa, Poland
 tel.: +48 34 39 06 735
 e-mail: bhp@semex.pl
 www.hfsafety.pl

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