nitriLAB

CAT III

CE2777

EN ISO 374-5:2016

VIRUS

EN 421:2010

EN ISO 374-1:2016+A1:2018/TYPE C

EN ISO 21420:2020

Rękawice ochronne kat. III

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Rozmiary: XS/5-5.5, S/6-6.5, M/7-7.5, L/8-8.5, XL/9-9.5,

XXL/10-10.5, XXXL/11-11.5

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Opis produktu:

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**Sizes:** XS/5-5.5, S/6-6.5, M/7-7.5, L/8-8.5, XL/9-9.5,

XXL/10-10.5, XXXL/11-11.5

Packing: 200 pieces/pack, 10 pack/carton**Product description:**

Disposable nitrile gloves

Ambidextrous
Dexterity level 5**General**

Powder free nitrile examination gloves. The product is classed as Category III of Personal Protective Equipment (PPE) by the European PPE regulation EU 2016/425 and have been shown to comply with this Regulation through the Harmonized European Standard(s) such as EN ISO 21420:2020, EN ISO 374-1:2016+A1:2018, EN ISO 374-2:2019, EN 16523-1:2015+A1:2018 and EN ISO 374-5:2016.

Information provided does not reflect the actual duration of protection in the workplace due to other factors influencing the performance, such as temperature, abrasion, and degradation. Information regarding protection refers to the working surface, i.e. the 'palm' of the glove, which has been submitted to testing. Before use, gloves should be thoroughly checked for damage (cuts and holes). If any damage is found, avoid using the product. The information does not reflect the actual duration of protection in the workplace or the difference between mixtures and pure chemicals. This product does not protect against mechanical hazards. If the conditions of use require protection against radioactive contamination in accordance with EN 421:2010, ensure that the glove protrudes sufficiently beyond the sleeve of the garment. Chemical resistance has been assessed under laboratory conditions based on samples taken exclusively from the hand (except in cases where the glove is 400 mm or longer, in which case the cuff is also tested) and refers exclusively to the chemical substance tested. It may vary if the chemical is used in a mixture. It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation. When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use me significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves. Before use, inspect gloves for any defects or imperfections. The penetration resistance has been assessed under laboratory conditions and relates only to the tested specimen. Please follow the correct procedures when donning, doffing, adjusting the gloves preserving comfort and hygiene of the hand to avoid any contamination. Single use only.

Donning

Insert the glove to the hand carefully without damaging the glove.

Doffing

Grasp the outside of the glove from the wrist area. Peel the glove away from the hand, hold it in the opposite gloved hand. Slide an un-gloved finger under the wrist of the remaining glove, being careful not to touch the contaminated surface of the glove. Additional information will be supplied upon request.

Storage

Keep out from direct sunlight; store in a cool and dry place. To be stored in cool dry place away from sunlight. Storage procedures are the main factor in determining glove shelf life. Gloves should be kept in their packaging protected from sunlight, artificial light, and humidity and stored at temperatures between 5 °C -35 °C. Storage under these conditions should provide shelf life of three years.

Intended Purpose

A patient examination glove used for diagnostic and therapeutic procedures. Reduce risk of infection & provide general protection of both patient and user against germs, bacteria, blood and body fluids as well as prevention of cross contamination. The glove is intended to be used in Medical sector, Healthcare workers & medical staffs, Emergency medical service.

Disposal of packaging

All empty packagings should be disposed of in accordance with local, state, and federal regulations or recycled/reconditioned at an approved facility.

Warning: Compliance with the precautions and warnings in the instructions for use (IFU) of the medical devices ensures the safe use and application. This product contains accelerators and may cause allergic reactions in some individuals with Type IV allergies (Di-Thiocarbamate Type). If a reaction in the form of skin irritation is noticed, remove gloves immediately and wash the affected part with saline water. If there is no relief, seek medical advice.

The product is manufactured under a Quality Control System which has been satisfactorily assessed by SATRA Technology Europe Ltd. Bracetown Business Park, Clonee, Dublin, D15 YN2P, Ireland as meeting the requirement of PPE Regulation (EU) 2016/425 Module C2 and Type B certificate.

Protection against Chemical Hazards, according to EN ISO 374-1:2016+A1:2018/TYPE C

Chemical	Level	Mean Degradation(%)
40% Sodium Hydroxide (K)	Level 6	-30.6
30% Hydrogen Peroxide (P)	Level 1	-2.3
37% Formaldehyde (T)	Level 3	-8.4

EN ISO 374-1:2016+A1:2018 performance levels are as follows

Performance Level	0	1	2	3	4	5	6
Breakthrough time (min)	<10	>10	>30	>60	>120	>240	>480

EN ISO 374-5:2016

Protection against bacteria and Fungi Pas
Protection against virus Pas

Complies with standard EN 421:2010

Food contact

Approved for contact with food in accordance with EU Commission Regulation No. 10/2011. Suitable for contact with dry, acidic and alcoholic food products, not suitable for contact with fats and oils.

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